



## GEA 600 LITER HIGH SHEAR MIXER, MODEL PMA 600 ADVANCED

SKU: 2908-1

*Unused* 600 liter GEA high shear mixer, model PMA 600 Advanced, 316L stainless steel product contact surfaces, 240-450 liter working capacity, 36 kw main impeller, 1-170 rpm, 11 kw chopper, 1000-3000 rpm, jacketed bowl, machine designed to 16 bar shock rating, full vacuum rated bowl, 8" side discharge with valve, integrally mounted Quadro Comil, model U20, integrally mounted blade hoist, associated control systems with Allen Bradley plc and HMI designed for class 1 div 2 area, 480 volt, with GEA temperature control unit, model PMS 600, WIP design, through the wall design mounted on base with PL600 IBC hoist with (4) approximately 750 liter IBC's, PMA 600 Advanced serial# 14369.1310, built 2015 and unused.

### Product Specifications

*Chamber Volume: 600*

*Drawings: Yes*

*Manufacturer: GEA*

*Material Of Construction: Stainless Steel*

*Model: PMA 600 Advanced*


*Motor HP: 48*

*Serial Number: 14369.1310*

*Year Of Construction: 2015*

### PLEASE CONTACT:

Adam Covitt  
adam@fedequip.com  
fedequip.com  
216.271.3500





## All the machines and training you need under one roof



Federal Equipment Company buys and sells pharmaceutical manufacturing and processing equipment. The company's inventory includes equipment for the manufacture of pharmaceuticals, vitamins, nutraceuticals, supplements, ingredients, and food and beverage manufacturing. We generally stock equipment to produce solid dose, tablet and capsules, liquids, creams, gels, and complex aseptic processing and fill-finish operations, as well as solid dose (tablet and capsule) manufacturing.



Our extensive inventory is stored in more than 200,000 square feet of clean, controlled warehouse space which includes a training and formulation facility operated by Techceuticals.

### Solid Dose Processing / Manufacturing


Handling  
Granulation  
Tableting

- Dedusting
- Friability Testing
- Hardness Testing
- Metal Detecting
- Tablet Presses
- Printing
- Associated Turrets/Parts

Capsule Filling

- Banding
- Checkweighing
- Polishing
- Tooling/Change Parts
- Hard Gel Capsule Filling
- Soft Gel Capsule Filling

Coating  
Washing



### Aseptic Processing and Freeze Drying

Sterilizing  
Aseptic Processing  
Freeze Dryers / Lyophilizers

### All Packaging: Primary and Secondary

Techceuticals is a worldwide provider of solutions for tablet and capsule manufacturing operations; developed and tested over many years with proven results.

### SHOW SPECIAL:

Mention this ad when purchasing your next machine and receive one free pass to your choice of:

- The Manufacturing Process: Solid Dose Operations
- Tablet Press Operations Class